

Catheters for urodynamic



Products Catalogue



Catheters for urodynamic



Catheter for women 2-way

Women 2-way catheter is made of PVC and is available in four sizes with straight finish. The catheter is supplied sterile and included two connecting catheter. The set is intended for single use. The catheter was developed in cooperation with Czech and Slovak doctors.



Parameters:

- material PVC
- 2 side eyes
- straight finish
- connectors Luer Lock negative
- sterile
- 2-connecting hose 1200mm

Ordering data:

Code	Diameter [Fr.]	Length [mm]
A-Z2-06	6	400
A-Z2-08	8	300
A-Z2-12	12	300
A-Z2-16	16	300




Catheters for urodynamic



Catheter for women 3-way

Women 3-way catheter is made of PVC and is available in three sizes with straight finish. The catheter is supplied sterile and included three connecting catheter. The set is intended for single use. The catheter was developed in cooperation with Czech and Slovak doctors.



Parameters:

- material PVC
- 3 side eyes
- straight finish
- connectors Luer Lock negative
- sterile
- 3-connecting hose 1200mm

Ordering data:

Code	Diameter [Fr.]	Length [mm]
A-Z3-08	8	300
A-Z3-12	12	300
A-Z3-16	16	300



Catheters for urodynamic



Catheter for men 2-way

Men 2-way catheter is made of PVC and is available in three sizes with tiemann finish. The catheter is supplied sterile and included two connecting catheter. The set is intended for single use. The catheter was developed in cooperation with Czech and Slovak doctors.



Parameters:

- material PVC
- 2 side eyes
- tiemann finish
- connectors Luer Lock negative
- sterile
- 2-connecting hose 1200mm

Ordering data:

Code	Diameter [Fr.]	Length [mm]
A-M2-06	6	400
A-M2-08	8	360
A-M2-12	12	360




Catheters for urodynamic



Catheter for men 3-way

Men 3-way catheter is made of PVC and is available in two sizes with tiemann finish. The catheter is supplied sterile and included three connecting catheter. The set is intended for single use. The catheter was developed in cooperation with Czech and Slovak doctors.



Parameters:

- material PVC
- 3 side eyes
- tiemann finish
- connectors Luer Lock negative
- sterile
- 3-connecting hose 1200mm

Ordering data:

Code	Diameter [Fr.]	Length [mm]
A-Z3-08	8	360
A-Z3-12	12	360




Catheters for urodynamic



Rectal catheter with balloon

Rectal catheter is made of PVC and is available in two sizes with straight finish. The catheter is supplied sterile and included one connecting catheter. The set is intended for single use. The catheter was developed in cooperation with Czech and Slovak doctors.



Parameters:

- material PVC
- 2 side eyes
- straight finish
- connectors Luer Lock negative
- sterile
- 1-connecting hose 1200mm
- silicone balloon

Ordering data:

Code	Diameter [Fr.]	Length [mm]	Balloon [ml]
B-REK-08	8	380	20
B-REK-18	18	380	50

